Phonak Target Track.

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American Academy of Audiology Task Force on Pediatric Amplification (2013). Pediatric Amplification. American Academy of Audiology Clinical Practice Guidelines. <https://www.audiology.org/practiceguideline/clinical-practice-guidelines-pediatric-amplification/>

Bagatto, M., Scollie, S., Moodie, S. T., Seewald, R., Hyde, M., Glista, D., ... & Parsa, V. (2019). Protocol for the provision of amplification. Retrieved from <https://ir.lib.uwo.ca/nca/5/>

Cote, C. (2021). Applications of Data Analytics in Health Care. Business Insights. Retrieved from <https://online.hbs.edu/blog/post/data-analytics-in-healthcare>

Feilner, M., Rich, S., & Jones, C. (2016). Automatic and directional for kids - Scientific background and implementation of pediatric optimized automatic functions. Phonak Insight, retrieved from www.phonak.com/evidence

Flanigan, B., Lockwood, M., & Chang, C. (2017). Analytics to improve outcomes and reduce cost: Health systems and health plans can work together to jointly win the shift from volume to value. Deloitte. Retrieved from www.deloitte.com/us/population-health-analytics

Healthcare Data Analytics. (2023). Mercury Healthcare. Retrieved from <https://www.mercuryhealthcare.com/faq/what-is-healthcare-analytics-on-27/08/2024>

McCreery, R. W., Walker, E. A., Spratford, M., Bentler, R., Holte, L., Roush, P., ... & Moeller, M. P. (2015). Longitudinal Predictors of Aided Speech Audibility in Infants and Children. *Ear and Hearing*, 36, 245-375.

Phonak Pediatric Solutions. Product and age recommendation overview (2023). Retrieved from www.phonak.com/evidence/Age-Consideration-Guide

Tomblin, J. Bruce; Moeller, Mary Pat. Editorial: The Outcomes of Children with Hearing Loss Study. *Ear and Hearing* 36:p 1S-3S, November/December 2015. | DOI: 10.1097/AUD.0000000000000220

Visram, A. S., Roughley, A. J., Hudson, C. L., Purdy, S. C., & Munro, K. J. (2021). Longitudinal changes in hearing aid use and hearing aid management challenges in infants. *Ear and hearing*, 42(4), 961-972.

Wolfe, J., Duke, M., Schafer, E., Jones, C., & Rakita, L. (2017). Evaluation of Adaptive Noise Management Technologies for School-Age Children with Hearing Loss. *Journal of the American Academy of Audiology*, 28(5), 415-435.

World report on hearing. Geneva: World Health Organization; 2021. License: CC BY-NC-SA 3.0 IGO.

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