Phonak Insight

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Abrams, H. B. & Kihm, J. (2015). An introduction to MarkeTrak IX: A New Baseline for the Hearing Aid Market. *Hearing Review*, 22(6).

Appleton, J., & Voss S.C. (2020) Motion-based beamformer steering leads to better speech understanding and overall listening experience. *Phonak Field Study News*. Retrieved from: www.phonakpro.com/evidence, accessed August 23rd, 2022.

Appleton, J. (2020) AutoSense OS 4.0 - significantly less listening effort and preferred for speech intelligibility. *Phonak Field Study News* retrieved from www.phonakpro.com/evidence, accessed August 23rd, 2022.

Appleton, J. & König, G. (2014). Improvements in speech intelligibility and subjective benefit with binaural beamformer technology. *Hearing Review*, 21(11), 40-42.

Appleton, J. (2022). What Is Important to Your Hearing Aid Clients...and Are They Satisfied? *Hearing Review*. 29 (6).

Derleth, P., Georganti, E., Latzel, M., Courtois, G., Hofbauer, M., Raether, J., & Kuehnel, V. (2021). Binaural Signal Processing in Hearing Aids. *Seminars in Hearing*, 42, 206 - 223.

Economy, P. (2015). 26 Brilliant Quotes on the Super Power of Words. Retrieved from: <https://www.inc.com/peter-economy/26-brilliant-quotes-on-the-super-power-of-words.html>, accessed August 23rd, 2022.

Edwards, B. (2016). A Model of Auditory-Cognitive Processing and Relevance to Clinical Applicability. *Ear and Hearing*, 37(suppl.1), 85S-91S. Retrieved from: <https://doi.org/10.1097/AUD.0000000000000308>, accessed 23rd August 2022

Jensen N. S., Neher T., Laugesen S., Johannesson, R. B. & Kragelind, L. (2013). Laboratory and field study of the potential benefits of pinna cue-preserving hearing aids. *Trends in Amplification*. 17 (3/4). 171-188.

Jespersen, C. T., Kirkwood, B. C., Groth, J. (2021). Increasing the effectiveness of hearing aid directional microphones. *Seminars in Hearing*. 42: 224-236.

Hayes, D. (2019). Speech detection by direction. *Unitron White Paper*. Retrieved from: https://www.unitron.com/au/en_au/learn/speech-detection-by-direction0.html, accessed August 23rd, 2022.

Keidser, G., O'Brian, A., Hain, J., McLelland, M., & Yeend, I. (2009). The effect of frequency-dependent microphone directionality on horizontal localization performance in hearing-aid users. *International Journal of Audiology*, 48(11), 789-803.

Killion, M. C. (1997). The SIN report: Circuits haven't solved the hearing-in-noise problem. *Hearing Journal*, 50(10), 28-32.

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Latzel, M., Mejia, J., Young, T., & Hobi, S. (2022). ActiveVent™ Receiver provides benefit of open and closed acoustics for better speech understanding in noise and naturalness of own voice perception. *Phonak Field Study News*. Retrieved from www.phonakpro.com/evidence, accessed August 23rd, 2022.

Latzel, M., & Appleton-Huber, J. (2015). StereoZoom – Adaptive behaviour improves speech intelligibility, sound quality and suppression of noise. *Field Study News*. Retrieved from www.phonakpro.com/evidence, accessed August 23rd, 2022.

Lewis, S., Crandall, C., Valente, M., & Horn, J. (2004). Speech perception in noise directional microphones versus frequency modulation (FM) systems. *Journal of the American Academy of Audiology*, 15, 426-439.

Picou, E., Aspell, E., & Ricketts, T. (2014). Potential benefits and limitations of three types of directional processing in hearing aids. *Ear and Hearing*, 35(3), 339-52.

Raether, J. (2005). Real Ear Sound - A simulation of the pinna effect optimizes sound localization also with open fittings. *Phonak Field Study News*. Retrieved from www.phonakpro.com/evidence, accessed August 23rd, 2022.

Ricketts, T. & Henry, P. (2002). Evaluation of an adaptive, directional-microphone hearing aid. *International Journal of Audiology*, 41:2, 100-112.

Schulte, M., Meis, M., Krüger, M., Latzel, M., Appleton-Huber, J. (2018). Significant increase in the amount of social interaction when using StereoZoom. *Phonak Field Study News*. Retrieved from: www.phonakpro.com/evidence, accessed August 23rd, 2022.

Shashkevich, A. (2019). The power of language: How words shape people, culture. *Stanford University Communications*. Stanford University. Retrieved from: <https://news.stanford.edu/2019/08/22/the-power-of-language-how-words-shape-people-culture/> accessed August 23rd, 2022.

Stewart, E., Rakita, L. & Drexler, J. (2019). *Phonak Compendium: StereoZoom Part 1: The benefit of wirelessly connected narrow directionality in Phonak hearing aids for speech intelligibility.*

Thibodeau, L. (2014). Comparison of speech recognition with adaptive digital and FM wireless technology by listeners who use hearing aids. *American Journal of Audiology*, 23(2), 201-210. Thibodeau L. M. (2020). Benefits in Speech Recognition in Noise with Remote Wireless Microphones in Group Settings. *Journal of the American Academy of Audiology*, 31(6), 404-411.

Vas, V., Akeroyd, M. A., & Hall, D. A. (2017). A data-driven synthesis of research evidence for domains of hearing loss, as reported by adults with hearing loss and their communication partners. *Trends in Hearing*, 21, 1-25.

Wagener, K., Vormann, M., Latzel, M., & Müller, H. (2018). Effect of hearing aid directionality and remote microphone on speech intelligibility in complex listening situations. *Trends in Hearing*, 22, 1-12.

Walden B. E., Surr R. K., Cord M.T., & Dyrland O. (2004). Predicting hearing aid microphone preference in everyday listening. *J Am Acad Audiol*. 15(5):365-96.

Winkler, A., Latzel, M., & Holube, I. (2016). Open Versus Closed Hearing-Aid Fittings: A Literature Review of Both Fitting Approaches. *Trends in hearing*, 20, 1-13.

Winneke, A., Latzel, M., Appleton-Huber, J. (2018). Less listening- and memory effort in noisy situations with StereoZoom. *Phonak Field Study News*. Retrieved from: www.phonakpro.com/evidence, accessed August 23rd, 2022.

Winneke, A., Schulte, M., Vormann, M. & Latzel, M. (2020). Effect of directional microphone technology in hearing aids on neural correlates of listening and memory effort: an electroencephalographic study. *Trends in Hearing*. 24:1-16. Retrieved from: <https://journals.sagepub.com/doi/pdf/10.1177/2331216520948410>, accessed August 23rd, 2022.

Xu, J., & Han, W. (2014). Improvement of Adult BTE Hearing Aid Wearers' Front/Back Localization Performance Using Digital Pinna-Cue Preserving Technologies: An Evidence- Based Review. *Korean Journal of Audiology*, 18, 97-104. Retrieved from: <https://pdfs.semanticscholar.org/5fc8/4f3c69b924cffe2c2b1d4e14aa52b6cbb0dc.pdf>, accessed August 23rd, 2022

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