


# Custom Earmold & Swim Plug Order Form


## Prelocked Tubing (TRS)

Standard for all soft earmolds.



## Earmold Areas

Diagram showing ear areas with names.



## Earmold Colours



Clear/Transparent (21)



Rose Tint (T)



Black (06)



Blue (07)



Purple (08)



Red (10)



Orange (11)



Green (17)



Yellow (20)



Beige (26)

## Swim Plug Colours



Blue (07)



Red (10)



Green (17)



Yellow (20)



Beige (26)



Light Blue (34)

# Custom Earmold & Swim Plug Order Form

Please note that changes cannot be made without calling.

Date: \_\_\_\_ / \_\_\_\_ / \_\_\_\_  
DD MM Year

PO#: \_\_\_\_\_ Attention to \_\_\_\_\_  
(at Phonak): \_\_\_\_\_

## STEP 1: General Information

Ship to Account Number:

Phone#:

Company Name:

Address:

Contact Name (Audiologist/Dispenser):

First Name:

Last Name:

Age of patient: \_\_\_\_\_  
(Pediatric orders only)

SHIP TO

Claim Type:  AADL  BCEHP  CNESST  Eastern Health  
 Greenshield  ISC (NIHB)  Manitoba Health  Nisga'a  
 ODSP  REGIE  Supp.Health  VAC, DND, RCMP  
 WCB \_\_\_\_\_ province  WSIB  
 Other: \_\_\_\_\_

Claim # (Required): \_\_\_\_\_  
(Please Specify)

L        
R        
250 500 1K 2K 3K 4K

Please include audiogram with all orders. Audiogram required for Acoustically Optimized Venting (AOV).

## STEP 2: Order Information

A) Please send earmold(s) only

B) Please send earmold(s) with BTE hearing instrument(s):

L Model: \_\_\_\_\_ Colour: \_\_\_\_\_

R Model: \_\_\_\_\_ Colour: \_\_\_\_\_

Please also send me the following accessories: \_\_\_\_\_

C) Swim Plugs (Please go to step 5)

## STEP 3: Configuration

Venting

AOV Acoustically Optimized Vent (Audiogram required)

OR

Vent Size (Customer specified):

L: \_\_\_\_\_mm  R: \_\_\_\_\_mm

Ear Texture

Soft  Avg.  Firm

Canal Length

Short Med. Long  
 L  R  L  R  L  R

## STEP 4: Earmold Style



Canal Mold

L  R



Canal Lock

L  R



Helix Lock

L  R



Semi Skeleton

L  R



Skeleton

L  R



Carved Half Shell

L  R



Carved Full Shell

L  R



Full Shell

L  R

## STEP 5: Materials and Colours

S = Standard

<input type="checkbox"/> Please indicate if Swim Plug is for a Lyric wearer		Colour										Lacquer		
		Clear (21)	Rose Tint (T)	Black (06)	Orange (11)	Purple (08)	Red (10)	Blue (07)	Green (17)	Yellow (20)	Beige (26)	Light Blue (34)	Glossy Finish	No Lacquer (Clear only)
Acrylic Hard (Lucite)	<input checked="" type="checkbox"/> L <input checked="" type="checkbox"/> R	S											S	
Silicone Soft 60	<input type="checkbox"/> L <input type="checkbox"/> R	S											S	
Silicone Very Soft 40	<input type="checkbox"/> L <input type="checkbox"/> R	S											S	
Swim Plug - Floatable Silicone	<input type="checkbox"/> L <input type="checkbox"/> R													

## STEP 6: Tube/Connector Details

- Soft earmolds come standard with prelocked tubing (TRS).
- Hard earmolds come standard with glued tubing.
- Tube/Connectors are only available in clear/transparent.
- 13M is standard for acrylic and 13T is standard for silicone.

Tubing		
13 Medium Tube 13M*	<input type="checkbox"/> L	<input type="checkbox"/> R
13 Thick Wall Tube 13T	<input type="checkbox"/> L	<input type="checkbox"/> R
13 Dry Wall Tube 13D*	<input type="checkbox"/> L	<input type="checkbox"/> R
13 Dry Thick Wall Tube 13W	<input type="checkbox"/> L	<input type="checkbox"/> R

\*Not available with Prelocked Tubing

SlimTube	Please indicate Tube length (00*/0/1/2/3)	
SlimTube 3.0 (Bolero B, Sky B-M, -P, Et -PR)	L	R
SlimTube 4.0 (Lumity, Paradise, Marvel BTEs)	L	R
Power SlimTube II (Sky B-SP Et -UP, Naida B-SP Et -UP)	L	R

\*\*00 only available in SlimTube 4.0

## Special Instructions:

\_\_\_\_\_

\_\_\_\_\_