

Phonak Ladegeräte – Kompatibilitätsübersicht.

Infinio · Lumity

Phonak
ChargerGo
RICSPH I



075-0057-12

Phonak
ChargerGo RIC I



075-0055-12

Phonak
Charger RIC I



075-0054-12

Phonak
Charger Virto I



075-0075-12

Phonak
Charger Ease



075-0080-10

Phonak
Slim Charger



075-3024-12

Phonak Charger
Case Go™



075-3017-12

Phonak Life
Charger™



075-3023-12


Phonak Charger
Combi BTE 2



075-0065-12

Audéo I-Sphere	•							
Audéo I-R		•	•					
CROS I-R		•	•					
Virto I-R				•				
Audéo L-R/RT					•			
Audéo L-RL						•	•	
Audéo L-R Fit					•			
CROS L-R					•			
Naída L-PR								•
Sky L-PR								•
Slim L-R						•		
Roger Focus II								•

Die Charger-Kompatibilität älterer Geräte finden Sie hier: www.phonak.de/charger

 Kompatibel mit Phonak Power Pack (071-0003)