

Phonak Charger – Kompatibilitätsübersicht.

Lumity · Paradise · Marvel

Phonak Charger Ease	Phonak Slim Charger	Phonak Charger Case Go™	Phonak Life Charger™	Phonak Charge and Care	Phonak Charger Case Combi 2	Phonak Charger Combi BTE 2
						
075-0080-10	075-3024-12	075-3017-12	075-3023-12	075-0041-12	075-0063-12	075-0065-12

	Kompaktes Ladegerät	Kompaktes Ladegerät	Integrierter Akku Induktives Laden	Induktives Laden	Elektrische Trocknung	Trocknung mit Trockenkapsel
Phonak Audéo L-R/RT	●					
Phonak Audéo L-RL (Life)			●	●		
Phonak Audéo L-R Fit	●					
Phonak CROS L-R	●					
Phonak Naída L-PR						●
Phonak Sky L-PR						●
Phonak Slim L-R		●				
Phonak Audéo P-R/RT					●	●
Phonak Audéo P-R Fit						●
Phonak Audéo P-RL (Life)			●	●		
Phonak CROS P-R					●	●
Phonak Naída P-PR						●
Phonak Audéo M-R/RT					●	●
Phonak Bolero M-PR						●
Phonak Sky M-PR						●
Roger Focus II						●

● Kompatibel mit Phonak Power Pack (071-0003)