

# Allegato CeruStop

Dati tecnici

Phonak Slim R

Phonak Audéo Lumity (R, RL, RT)

Phonak Audéo Paradise (R, RL, RT, 13T, 312)

Phonak Audéo Marvel (R, RT, 13T, 312, 312T)

Ricevitore P

Dati accoppiatore da 2 cm<sup>3</sup>

ANSI / ASA S3.22-2014 (R2020)  
IEC 60118-0:2022

Guadagno acustico

Range di frequenza	<100 Hz - 6300 Hz (CeruShield)
	<100 Hz - 5900 Hz (CeruStop)



Sonova AG · Laubisrütistrasse 28  
CH-8712 Stäfa · Switzerland  
www.phonak.com

A Sonova brand

**PHONAK**  
life is on