

Bijlage CeruStop

Technische gegevens

Phonak Slim R

Phonak Audéo Lumity (R, RL, RT)

Phonak Audéo Paradise (R, RL, RT, 13T, 312)

Phonak Audéo Marvel (R, RT, 13T, 312, 312T)

P Receiver

2 cm³ couplergegevens

ANSI / ASA S3.22-2014 (R2020)
IEC 60118-0 : 2022

Akoestische versterking

Frequentiebereik	<100 Hz - 6300 Hz (CeruShield)
	<100 Hz - 5900 Hz (CeruStop)



0459



Sonova AG · Laubisrütistrasse 28
CH-8712 Stäfa · Switzerland
www.phonak.com

A Sonova brand

PHONAK
life is on