

# Anhang CeruStop

Technische Daten

Phonak Slim R

Phonak Audéo Lumity (R, RL, RT)

Phonak Audéo Paradise (R, RL, RT, 13T, 312)

Phonak Audéo Marvel (R, RT, 13T, 312, 312T)

P Receiver

2-cm<sup>3</sup>-Kuppler-Daten

ANSI / ASA S3.22-2014 (R2020)  
IEC 60118-0:2022

---

## Akustische Verstärkung

|                 |                                |
|-----------------|--------------------------------|
| _____           | _____                          |
| Frequenzbereich | <100 Hz - 6300 Hz (CeruShield) |
| _____           | _____                          |
|                 | <100 Hz - 6000 Hz (CeruStop)   |

---



0459



Sonova AG · Laubisrütistrasse 28  
CH-8712 Stäfa · Switzerland  
[www.phonak.com](http://www.phonak.com)

A Sonova brand

**PHONAK**  
life is on