

Phonak ChargerGo RIC 1

Technical data.

| | |
|-------------------------------------------------------------|---------------------------------------------------------|
| Dimensions (L x W x H) | 79 x 40 x 55 mm (3.1 x 1.5 x 2.2") +/- 5% |
| Weight | 68 g (2.4 oz) +/- 10% |
| Color | Champagne |
| Operating conditions | +5° to +40° Celsius (+41° to +104° Fahrenheit) |
| Transport and storage conditions | -20° to +60° Celsius (-4° to +140° Fahrenheit) |
| Humidity during transportation and storage | 0% to 93% (non condensing) |
| Humidity during operation | 0% to 90% (non condensing) |
| Atmospheric pressure during use, storage and transportation | 500 hPa to 1060 hPa |
| Included Parts | Power supply USB-C cable |
| External Interfaces | USB-C port |
| Battery | Li-Ion (built-in) |
| Charging time | See User Guide |
| Input | Type: DC, stabilized 5 V +/- 5% Max. current: 350 mA |



The battery cannot be over-charged. It may become warm during charging.