

# Phonak Charger RIC 1

## Technical data.

Dimensions (L x W x H)	79 x 40 x 55 mm (3.1 x 1.5 x 2.2") +/- 5%
Weight	60 g (2.1 oz) +/- 10%
Color	Champagne
Operating conditions	+5° to +40° Celsius (+41° to +104° Fahrenheit)
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit)
Humidity during transportation and storage	0% to 93% (non condensing)
Humidity during operation	0% to 90% (non condensing)
Atmospheric pressure during use, storage and transportation	500 hPa to 1060 hPa
Included Parts	Power supply USB-C cable
External Interfaces	USB-C port
Charging time	See User Guide
Input	Type: DC, stabilized 5 V +/- 5% Max. current: 350 mA

