

Phonak cShell 6.0 order form.

Customer information

Ship to account number: _____

Address: _____

City: _____ State: _____ Zip: _____

Bill to account number: _____

Third party patient number: _____

Fitting date: _____

Purchase order number: _____

Contact name: _____

Phone: _____

Email: _____

Hearing instrument

Qty	Instrument type	Technology level			
		90	70	50	30
	Phonak Audéo™ I Sphere	<input type="checkbox"/>	<input type="checkbox"/>	-	-
	Phonak Audéo I-R	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Phonak CROS I-R	-	-	-	-

Specify color choice: _____

See page 2 for color options

cShell style



cShell Acrylic (062-0056-01)

L R

Shell colors: Black (06), Cocoa (22), Pink (26), Tan (14), Brown (28), Red (R), Blue (B), Transparent (21)

Faceplate colors: Black (06), Cocoa (22), Pink (26), Tan (14), Brown (28)

Accessories and Roger™

	TV Connector (076-3006-0611)	<input type="checkbox"/>
Accessories	PartnerMic (076-4001-P511)	<input type="checkbox"/>
	RemoteControl (076-0065-P5)	<input type="checkbox"/>
Roger	Roger On™ 3 (056-3032-xx011) Champagne (P5), Graphite Gray (P7)	<input type="checkbox"/> P5 <input type="checkbox"/> P7
	Roger Select™ 3 (056-3033-xx011) Champagne (P5), Graphite Gray (P7)	<input type="checkbox"/> P5 <input type="checkbox"/> P7
	Roger Table Mic (056-4012-P5011)	<input type="checkbox"/>

Preferences

If necessary, may we change the following: Please call

Yes	No	
<input type="checkbox"/>	<input type="checkbox"/>	UP to P OR P to M OR M to S if audiogram permits
<input type="checkbox"/>	<input type="checkbox"/>	Build larger if required
<input type="checkbox"/>	<input type="checkbox"/>	I would like to receive a paper user guide

Special Instructions

Patient information

Last name: _____

First name: _____

Age: _____ Gender: _____

Is this patient order being fully or partially paid for by state or federal government funds? Yes No

Audiogram (Required for AOV):

	Hz	250	500	1K	2K	3K	4K
Left:	AC						
Right:	AC						

HI warranty

4th year 5th year (Patient age required at time of order)

Shipping and handling options[†]

2 Business Days (\$20.99) C.O.D. (\$21.99 + shipping)

3-5 Business Days (\$19.99)

[†] Based on location and courier

Options

Shell Color	Standard color (Transparent 21)	<input checked="" type="checkbox"/> S
	Custom color	Code: _____
Faceplate Color	Standard color	<input checked="" type="checkbox"/> S
	Custom color	Code: _____
Receiver	S Receiver	L <input type="checkbox"/> R <input type="checkbox"/>
	M Receiver	L <input type="checkbox"/> R <input type="checkbox"/>
	P Receiver	L <input type="checkbox"/> R <input type="checkbox"/>
	UP Receiver (may result in cShell XL)	L <input type="checkbox"/> R <input type="checkbox"/>
Length	0/1/2/3	Left: _____ Right: _____
	AOV (enter audiogram)	<input checked="" type="checkbox"/> S
Vent size	Customer specific 12/20/25/3L/none	Left: _____ mm Right: _____ mm
	Standard wax system**	<input checked="" type="checkbox"/> S
Wax system	Custom wax system	
	Helix Lock (HL)	<input type="checkbox"/>
Retention*	Canal Lock (CL)	<input type="checkbox"/>
	Skeleton Lock (SL)	<input type="checkbox"/>
	Transparent (21)	<input type="checkbox"/>
Removal line (RF)	Line extension + customer specific	
	Extra retention (ERS)	<input type="checkbox"/>

* Not available for cShell Titanium

** HF4 for UP receiver, CeruStop for all other receivers

Additional Information

Alternative cShell Build Style

We offer an alternative cShell build style with an anterior tube outlet. This build style provides high cosmetic acceptance in spite of a slightly larger size and can be used in case of build problems due to anatomical limitations. This build style is available for acrylic shells only.



Fig.1 cShell - standard tube




Fig.2 cShell - anterior tube outlet

Acoustically Optimized Vent (AOV)

We recommend to select AOV. We will set the correct vent diameter and style consider audiogram/shell geometry/feedback canceler and required low frequency gain. The client's audiogram is mandatory for the calculation of the AOV.

Extra Retention Surface (ERS)

We offer a rougher surface for cShell titanium. The extra retention surface can improve retention in the ear.

Wax systems

Sensor cShell	cShell Acrylic	Description
X	X	No wax guard
CS	CSB	CeruStop
H4	H4B	HF4 Filter, significantly smaller than HF3

Hearing instrument compatibility


The Sensor cShell 6.0 is exclusively provided for Audéo Fit. The cShell 6.0 acrylic is compatible to Infinio hearing aids or newer. It is not backward compatible to Lumity or earlier generation hearing aids.

Color codes

Code	Color
06	Black
13/21	Transparent
14	Tan
22	Cocoa
26	Pink
28	Brown
R/B	Red/Blue
U0	Titanium Gray

Receiver tube length

Measure the tube length to select the correct receiver. Measure both ears for binaural fitting.



Infinio hearing instrument color options

