

Phonak SlimTip/CROS Tip 6.0 order form.

Customer information

Ship to account number: _____

Address: _____

City: _____ State: _____ Zip: _____

Bill to account number: _____

Third party patient number: _____

Fitting date: _____

Purchase order number: _____

Contact name: _____

Phone: _____

Email: _____

Patient information

Last name: _____

First name: _____

Age: _____ Gender: _____

Is this patient order being fully or partially paid for by state or federal government funds? Yes No


Hearing instrument

| Qty | Instrument type | Technology level | | | |
|-----|------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| | | 90 | 70 | 50 | 30 |
| | Phonak Audéo™ I Sphere | <input type="checkbox"/> | <input type="checkbox"/> | - | - |
| | Phonak Audéo I-R | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | Phonak CROS I-R | - | - | - | - |

Specify color choice: _____

See page 2 for color options

SlimTip style




SlimTip Acrylic (062-0063-01)

L R

Shell style: Hollow

Standard color: Transparent (21)

Custom colors: Tan (14), Cocoa (22), Pink (26), Brown (28), Red/Blue (R/B)




SlimTip Silicone (062-0064-01)

L R

Shell style: Solid

Standard color: Transparent (21)




SlimTip Titanium AV (062-0065-01)

Only available with Active Vent Receiver MAV

L R

Shell style: Hollow

Color: Titanium gray (U0)



CROS Tip (062-0059)

L R

Shell style: Hollow

Standard color: Transparent (21)

Custom colors: Tan (14), Cocoa (22), Pink (26), Brown (28), Red/Blue (R/B)

Preferences

If necessary, may we change the following: Please call

Yes No

 P to M **OR** M to S if audiogram permits

 Build larger if required

 I would like to receive a paper user guide

Audiogram (Required for AOV):

| | Hz | 250 | 500 | 1K | 2K | 3K | 4K |
|--------|----|-----|-----|----|----|----|----|
| Left: | AC | | | | | | |
| Right: | AC | | | | | | |

HI warranty

4th year 5th year (Patient age required at time of order)

Shipping and handling options[†]

2 Business Days (\$20.99) C.O.D. (\$21.99 + shipping)

3-5 Business Days (\$19.99)

[†] Based on location and courier

Options

| | | |
|-------------------------|---|---|
| Shell Color | Standard color | <input checked="" type="checkbox"/> S |
| | Custom color | Code: _____ |
| Connection | S Receiver | <input type="checkbox"/> L <input type="checkbox"/> R |
| | M Receiver | <input type="checkbox"/> L <input type="checkbox"/> R |
| | P Receiver | <input type="checkbox"/> L <input type="checkbox"/> R |
| Length | Length 00*/0/1/2/3 | Left: _____ Right: _____ |
| | AOV (enter audiogram)** | <input checked="" type="checkbox"/> S |
| Vent size | Customer specific** | Left: _____ mm Right: _____ mm |
| | Cavity (Y) | <input checked="" type="checkbox"/> S |
| Vent style | Regular (P) | <input type="checkbox"/> |
| | Helix Lock*** (HL) | <input type="checkbox"/> |
| Retention** | Canal Lock (CL) | <input type="checkbox"/> |
| | Skeleton Lock (SL) | <input type="checkbox"/> |
| | Removal line (always included) | <input checked="" type="checkbox"/> |
| Removal line | Line extension + customer specific length | |
| | Extra retention (ERS) | <input type="checkbox"/> |
| Titanium surface | TV Connector | <input type="checkbox"/> |
| | PartnerMic | <input type="checkbox"/> |
| | Remote Control | <input type="checkbox"/> |
| Accessories | Roger On™ 3 (056-3032-xx011) Champagne (P5), Graphite Gray (P7) | <input type="checkbox"/> P5 <input type="checkbox"/> P7 |
| | Roger Select™ 3 (056-3033-xx011) Champagne (P5), Graphite Gray (P7) | <input type="checkbox"/> P5 <input type="checkbox"/> P7 |
| | Roger Table Mic 3 (056-4012-P5011) | <input type="checkbox"/> |

Special Instructions


Additional Information

Acoustically Optimized Vent (AOV)

We recommend to select AOV. We will set the correct vent diameter and style consider audiogram/shell geometry/feedback canceler and required low frequency gain. The patient's audiogram is mandatory for the calculation of the AOV.


Vent styles

We offer two vent options for SlimTips.



Regular:

- Classic vent solution
- Fixed diameter (not enlargeable)




Cavity:

- Optimized for best fitrate
- Easy to enlarge by the fitter

Tube length

Measure the tube length to select the correct SlimTube. Measure both ears for binaural fitting.



Wax Systems

The following wax systems are used for different types of SlimTip/CROS Tip.

| | |
|---------------------|------------------------------------|
| SlimTip Acrylic | CeruStop mounted in SlimTip* |
| SlimTip Silicone | CeruStop in receiver |
| SlimTip Titanium AV | HF3 mounted in ActiveVent receiver |
| CROS Tip | No wax filter needed |

* CeruStop from receiver must be removed

SlimTip Titanium AV

This variant of the SlimTip is designed for ActiveVent™ (AV) and is compatible to the MAV receiver (MAV = Medium ActiveVent) only. There is no additional vent option needed. It is offered in durable Titanium to maximize fitrate and increase comfort.

Our SlimTip Link AV connector makes the assembly of the MAV Receiver simple and easy.

ActiveVent is compatible to Infinio rechargeable RIC hearing aids.

Extra Retention Surface (ERS)

We offer a rougher surface for Titanium. The extra retention surface can improve retention in the ear.

Hearing instrument compatibility

The SlimTip 6.0 portfolio is compatible to Infinio hearing aids or newer. It is not backward compatible to Lumity or earlier generation hearing aids.

Color codes

| Code | Color |
|------|---------------|
| 14 | Tan |
| 21 | Transparent |
| 22 | Cocoa |
| 26 | Pink |
| 28 | Brown |
| R/B | Red/Blue |
| U0 | Titanium Gray |

Infinio hearing instrument color options

