

# Anexo CeruStop

Dados técnicos

Phonak Slim R

Phonak Audéo Lumity (R, RL, RT)

Phonak Audéo Paradise (R, RL, RT, 13T, 312)

Phonak Audéo Marvel (R, RT, 13T, 312, 312T)

Receptor P

Dados do acoplador de 2 cm<sup>3</sup>

ANSI/ASA S3.22-2014 (R2020)  
IEC 60118-0: 2022

## Ganho acústico

Faixa de frequência	<100 Hz-6300 Hz (CeruShield)
	<100 Hz - 6000 Hz (CeruStop)



0459



Sonova AG · Laubisrütistrasse 28  
CH-8712 Stäfa · Switzerland  
www.phonak.com

A Sonova brand

**PHONAK**  
life is on