

# Phonak CROS P Order Form

Phonak CROS P is a solution for unilateral hearing loss designed to help the wearer follow conversations. CROS P is worn on the unaidable ear to pick up sounds and wirelessly transmits them to a Phonak Audéo™ P hearing aid worn on the ear with better hearing. CROS P has no acoustical output.

## Step 1: Order details

Account Number: 20  
 Hearing Clinic:  
 Delivery Address:  
 Contact Name: Phone Number:  
 Email Address:  
 Date Ordered: Date Required:

## Step 2: Patient information

Name:

## Step 3: VIP service

VIP 24 hour service (\$57.75 incl GST)

## Step 4: Choose hearing instrument

(see following page for HSP codes)

Model	Part No.	HSP	Colour	Qty
CROS P-13	050-0861-xx	N0275CR		


## Colour chart



## Step 5: Choose your Acoustic Coupling Option


### CROS Slim Tube 4.0 with Domes

CROS Slim Tube 4.0 is available in left/right and four different lengths (0, 1, 2, 3) with standard dome available in three sizes (S, M, L)




### CROS Slim Tube 5.0 with Slim Tip<sup>^</sup>

CROS Slim Tube 5.0 is available in left/right and four different lengths (0, 1, 2, 3), with a custom made earpiece



### CROS Tip 4.0<sup>^</sup>

Fully custom made earpiece for CROS, available in four different cable lengths (0, 1, 2, 3)



L R Select options in section a and c only below

L R Select options in section c and d only below

L R Select options in section b and d only below

<sup>^</sup> Impression required

<sup>^</sup> Impression required

## Step 6: Choose from the following options

### (a) Domes 4.0

Dome style	Dome size	
	Left	Right
Cap	Left	Right
Open	Left	Right
Vented	Left	Right
Power	Left	Right

### (b) CROS Tip Cables – must be paired with a CROS Tip below



0-L 0-R  
 1-L 1-R  
 2-L 2-R  
 3-L 3-R

### (c) CROS Slim Tube 5.0

054-1180	CROS Tube 5 0L	054-1184	CROS Tube 5 0R
054-1181	CROS Tube 5 1L	054-1185	CROS Tube 5 1R
054-1182	CROS Tube 5 2L	054-1186	CROS Tube 5 2R
054-1183	CROS Tube 5 3L	054-1187	CROS Tube 5 3R

### (d) Slim Tip 4.0 or CROS Tip

	Shell colour	Clear (21)	Rose Tint (1)	Beige/Pink (26)	Tan (14)	Vent style	Cavity	Regular	Shell attribute	Hollow	Solid
SlimTip Acrylic (062-0006-01)											
SlimTip Silicone (062-0007-01)											
CROS Tip 4.0 Hard Acrylic (062-0016-**) - must be paired with a CROS Tip cable											

### Generic options (all earpieces)

Removal line Transparent (21) Removal line extension + 5mm

Locks Canal Skeleton

Vent size AOV

L: \_\_\_mm R: \_\_\_mm

S = Standard, will default if no other selection ticked

# Phonak CROS P Order Form

Phonak CROS P is a solution for unilateral hearing loss designed to help the wearer follow conversations. CROS P is worn on the unaidable ear to pick up sounds and wirelessly transmits them to a Phonak Audéo™ P hearing aid worn on the ear with better hearing. CROS P has no acoustical output.

---

## Step 7: Choose your wireless Paradise hearing instrument for the better ear

Audéo P-13 (NOTE: Please also complete the Phonak Audéo P-RIC order form to specify correct inclusions).

---

## Step 8: Choose your accessories

### Wireless Communication Accessories

Please use the Phonak Audéo P-RIC order form to specify accessories.

### Roger Accessories

Refer to Phonak Roger Order Form

---

## Step 9: Check List

Phonak Audéo P RIC order form accompanies this order

Roger order form accompanies this order (if applicable)

---

## Additional information

### HSP code

Model	HSP codes	Category
CROS P-13	N0275CR	NS

---

## Models



CROS P-13