

Phonak Terra+™ Order Form

Phonak Terra+ provide a balanced set of features for a personalised hearing experience with RIC and BTE.



Step 1: Order details

Account Number: _____
 Hearing Clinic: _____
 Delivery Address: _____
 Contact Name: _____ Phone Number: _____
 Email Address: _____
 Date Ordered: _____ Date Required: _____

Step 2: Patient information

Name: _____
 Audiogram: (Required for AOV vent option)

	kHz	0.25	0.5	1	2	4
LEFT AC (dB HL)						
RIGHT AC (dB HL)						

Step 3: Warranty/VIP service

4th year warranty (charges will apply)

VIP 24 hour service (\$57.75 incl GST)

Step 4: Choose hearing instrument (see following page for HSP codes and colour options)

		RIC		
L	R	Model	Part No.	Colour
		Terra+ RIC-312	050-1107-xx	

		BTE		
L	R	Model	Part No.	Colour
		Terra+ BTE-M	050-1104-xx	
		Terra+ BTE-SP	050-1105-xx	
		Terra+ BTE-UP	050-1106-xx	

Step 5: Choose Sound Delivery System 5.0 for Terra+ RICs only

Standard		Medium		Power		UltraPower – cShell only HF3 filter recommended	
S 0-L	S 0-R	M 0-L	M 0-R	P 0-L	P 0-R	UP 0-L	UP 0-R
S 1-L	S 1-R	M 1-L	M 1-R	P 1-L	P 1-R	UP 1-L	UP 1-R
S 2-L	S 2-R	M 2-L	M 2-R	P 2-L	P 2-R	UP 2-L	UP 2-R
S 3-L	S 3-R	M 3-L	M 3-R	P 3-L	P 3-R	UP 3-L	UP 3-R

Dome style		
Cap	Left	Right
Open	Left	Right
Vented	Left	Right
Power	Left	Right

Dome size		
Small	Medium	Large
Small	Medium	Large
Small	Medium	Large



SimTip (removeable earpiece)	
L	R
Material	
Hard Acrylic	Hard Acrylic
Soft Silicone	Soft Silicone
Shell colour	
<input type="checkbox"/> Clear (21)	<input type="checkbox"/> Clear (21)
Rosa Tint (T)	Rosa Tint (T)
Beige/Pink (26) (hard only)	Beige/Pink (26) (hard only)
Tan (14) (hard only)	Tan (14) (hard only)

cShell Acrylic	cShell Titanium*^
L	R
Shell colour	
<input type="checkbox"/> Beige/Pink (26)	<input type="checkbox"/> Beige/Pink (26)
Clear (21)	Clear (21)
Tan (14)	Tan (14)
Faceplate colour	
<input type="checkbox"/> Pink (26)	<input type="checkbox"/> Pink (26)
Cocoa (22)	Cocoa (22)
Brown (28)	Brown (28)
Transparent (13)	Transparent (13)
Tan (14)	Tan (14)
Black Transp. (47)	Black Transp. (47)

Impressions

New impressions enclosed

Use scanned data
Existing serial number: _____

SlimTip

cShell Titanium

cShell

Okay to change to cShell if SlimTip not possible

* Please note cShell Titanium products are not included free of charge with hearing aid order.

S = Standard, will default if no other selection ticked

^ Please note cShell Titanium is only available in Titanium shell colour and Black faceplate colour.

Venting (mm)*

L	R
AOV	
Smart vent size calculation. Requires AC thresholds.	
Occluded	
1.0mm	
1.1 - 1.6mm	
1.7 - 2.2mm	
2.3 - 3.0mm	
Open	

*Cerustop default option for SlimTip HA and only option for SlimTip Titanium.

Please note: A coupling code representing actual vent size is printed on custom earpieces. Please enter this code into Phonak Target when programming your client's hearing aids.

Wax management

L	R	SlimTip Hard Acrylic
<input type="checkbox"/>	<input type="checkbox"/>	Cerustop (CS)
Ext. receiver tube + wax spring (WE)		
None (X)		
L	R	SlimTip Soft Silicone
<input type="checkbox"/>	<input type="checkbox"/>	None (X)
L	R	cShell
<input type="checkbox"/>	<input type="checkbox"/>	Cerustop (CS)
HF3 (recommended for UP Rx) (H3)		
HF4 (H4)		
Ext. receiver tube + wax spring** (WE)		

**Not applicable for cShell Titanium.

Handling

L	R
	Removal Line (RF)
	Removal Line + 5mm extension (RFE)
	Canal Lock (CL)
	Skeleton Lock (SL)

Step 6: Choose Sound Delivery System for Terra+ BTEs only

Tone Hook with Earmould
Tone Hooks are standard with order. If required, please fill out the Earmould 4.0 Order Form

SlimTube with domes or SlimTip
Please complete the section below

Terra+ BTE-UP is not compatible with SlimTube. Please fit with Earmould only



SlimTube 4.0 Size

(SlimTube 4.0 is not compatible with Terra+ BTE-UP)

00-L	00-R
0-L	0-R
1-L	1-R
2-L	2-R
3-L	3-R

SlimTip 4.0

SlimTip Acrylic (062-0006-01)

SlimTip Silicone (062-0007-01)

S = Standard, will default if no other selection ticked

Connection type	Shell colour	Vent style	Cavity	Regular
SlimTip	Clear (21)	Rose Tint (7)	Beige/Pink (26)	Tan (14)

Vent size: AOV L: ___ mm
R: ___ mm

Removal line: Canal Skeleton

Locks: Canal Skeleton

Step 7: Choose accessories

TV Connector (076-3006-06)

Remote Control (076-0065-P5)

PartnerMic (076-4001-P510)

Special requests:

Additional information

HSP codes

Model	HSP codes
Terra+ RIC-312	B824PHO
Terra+ BTE-M	B827PHO
Terra+ BTE-SP	B826PHO
Terra+ BTE-UP	B825PHO

Models



Housing colours

P1	P3	P6	H0
Sand Beige	Sandalwood	Silver Gray	Beige

Accessories

